







Attends Cover-Dri is a range of surface protection products designed for short, medium and long term use. Cover-Dri has a **unique hexagon quilt design** with absorbent material in each pocket. This helps to spread the fluid between the hexagons and prevents leakage/spillage. Each quilted pocket ensures that there is protection across the product, rather than pooling in one area, and allows the product to be thin and effective. Each pocket contains **super absorbents**, which helps to draw and lock fluids away from the patient, and gives outstanding protection and containment and helps to reduce the risk of leakage.

A pale blue non slip back sheet, ensures that the protective pad stays in place. Fewer pads need to be used as Cover-Dri stays in place.

Product Range

Attends Cover-Dri is available in a number of sizes and absorbencies to help you select the right product for the right patient.

Product	Features	Benefits	Uses*
Cover Dri Plus 40x60cm 	Non Slip Waterproof Backsheet Individual absorbent cushions Lock away superabsorbent core All Round Leakage barrier White top sheet Virgin Fluff pulp	Remains in place to deliver maximum protection Quickly absorbs wetness whilst drainage channels rapidly distribute the urine, and prevents wrinkling of the core. Improves skin dryness and leakage protection Prevents leakage onto bed/chair Allows assessment of colour of fluids/discharge Reduced risk of cross infection	<ul style="list-style-type: none"> • Bed, Wheelchair and Chair protection • Procedures for example - Irrigation of small wounds • Light urinary incontinence • Incubators
Cover Dri Plus 60x60cm Cover Dri Super 60x60cm 	As above	As above	<ul style="list-style-type: none"> • Bed, and Chair protection - ideal for use in Theatres • Light urinary incontinence • Antenatal and Post natal care
Cover Dri Plus 60x90cm Cover Dri Super 60x90cm 	As above	As above	<ul style="list-style-type: none"> • Bed, and Chair protection • Use along length of bed, in delivery suite • Light/moderate urinary incontinence
Cover Dri Plus 80x170cm 	As above plus two wings	As above plus can be tucked under mattress to keep securely in place. Can also be used for cleaning the patient, when removing the product.	<ul style="list-style-type: none"> • Bed, and Chair protection, can be used to cover beds, for clients lying on top, during the day • Light/moderate urinary incontinence

* This list is not exhaustive, but indicates the range of uses these pads are available for.



Cover Dri's unique technology makes it the ideal product for procedures and containment. Unlike other underpads you may have used, if you leave this pad in situ for a couple of minutes, **fluid** will be **locked away** into the core, allowing you to remove the pad, and fluid, in one go.

Cover Dri is a product designed for clinical uses, and where non sterile containment is required. The virgin pulp reduces the chance of cross infection.

Cover-Dri pads can be used to protect bed and chair fabrics/surfaces.

Cover Dri Plus (5) is available in four sizes, 40x60cm, which is ideal for chairs, 60x60cm and 60x90cm are designed for beds and larger surfaces. 80x170cm has tuckable flaps, for improved fixation.

CoverDri comes in easy to use bags, of 50 products (or 30 for the 80x170cm), also reducing the risk of contamination

Uses	Needs	Products
Barium Enema	Contain, oily liquid, and be changed regularly	Wadding Pad/Tuckable 80x170cm
Body Fluid Leakage	Absorption, containment, and white, to be able to see changes in fluid i.e. merconium staining	Full range of Cover Dri
Faecal Containment	Containment	Full range of Cover Dri
General Theatre Use	Absorption, containment of body fluids, virgin pulp, to prevent cross infection	Cover Dri 60x60cm, or 60x90cm
Irrigation of wounds	Absorption, containment	Cover Dri 40x60cm, 60x60cm, or 60x90cm
Light Urinary Incontinence	Absorption, containment, for patients with a low output/leakage	Full range of Cover Dri
Obstetrics & Gynaecology	Containment of body fluids or cleansing fluids	Cover Dri 60x60cm, 60x90cm,
Procedures	Absorption, containment	Full range of Cover Dri
Special Care Baby Unit	Absorption of urine, and containment of faeces, soft top sheet for comfort	Cover Dri 40x60cm

For more information visit www.attendshealthcare.co.in
 E-mail: attends@walnutgroup.net or Customer care: 1860 425 2211