

# Attends® Pull-Ons

**Question:**

**When does the bladder weakness occur?**

**Answer:**

The client has reduced bladder and/or bowel sphincter control and incontinence can occur without the client being aware.

**What happens?**

Loss of urine and sometimes also faeces, several times during the day and night.

Attends Pull-Ons is a highly absorbent pant that can be pulled up and down by the client or the caregiver. They can manage both urinary and faecal incontinence. The pants are white for discretion.

Due to their flexible shape with asymmetric leg-cut, Attends Pull-Ons give the feeling of normal underwear.

Attends Pull-Ons have a highly absorbent core and anti-leakage barriers, which help to prevent leakage. The urine is kept away from the skin and the odour protection core locks unpleasant smells inside the product.

A wide elastic panel in the waistband (both front and back) and soft leg elastics provide optimum comfort for the wearer. The cloth-like backsheet of Attends Pull-Ons is made of soft, breathable material, which contributes to healthy skin. If needed, the side panel can easily be torn open, for quick disposal.

Attends Pull-Ons are available in a range of sizes.



| Size                     |    | Small                        | Medium                       | Large                        | X-Large                      | X-Small                      |
|--------------------------|----|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|
| Waist size               | cm | 56-75                        | 75-110                       | 100-160                      | 120-175                      | 45-60                        |
| Weight                   | g  | 53                           | 60                           | 64                           | 81                           | 63                           |
| Abs. (ISO- 11948-1)      | ml | 993                          | 1079                         | 1095                         | 1412                         | 1381                         |
| Absorption Core Material |    | Fluff & Superabsorbent       | Fluff & Superabsorbent       | Fluff & Superabsorbent       | Fluff & Superabsorbent       | Fluff & Superabsorbent       |
| Topsheet Material        |    | Non Woven                    | Non Woven                    | Non Woven                    | Non Woven                    | Non Woven                    |
| Backsheet Material       |    | Breathable textile backsheet | Breathable textile backsheet | Breathable textile backsheet | Breathable textile backsheet | Breathable textile backsheet |
| Leg Seam Material        |    | Lycra                        | Lycra                        | Lycra                        | Lycra                        | Lycra                        |
| Waist Band Material      |    | Lycra                        | Lycra                        | Lycra                        | Lycra                        | Lycra                        |
| Odour Protection         |    | Yes                          | Yes                          | Yes                          | Yes                          | Yes                          |
| Pack Count               |    | 1 x 21                       | 1 x 14                       | 1 x 14                       | 1 x 14                       | 1 x 21                       |
| Case Count               |    | 4 x 21                       | 4 x 14                       | 4 x 14                       | 4 x 14                       | 4 x 21                       |
| Cases per pallet         |    | 30 x 4 x 21                  | 30 x 4 x 14                  | 30 x 4 x 14                  | 30 x 4 x 14                  | 24 x 4 x 21                  |